

FLUORESCENCE IMMUNOASSAY ANALYSER

Features

- Quantitative/Qualitative Results
- 4.3" LCD touch Screen
- Rechargeable Li-ion battery
- LIS, HIS compatibility



Test Parameters:

Troponin I, NT-proBNP, CK-MB, H-FABP, Myoglobin, CRP, Procalcitonin, D-dimer, TSH, PSA, AFP, CEA, mAlb, CysC, HbA1c, Strep A

Test Parameters

Catalogue No.	Product Name	Sample	Test Range	Kit size	CE status
FIA-AFP-001	AFP FIA Rapid Test Cassette	S/P	5~400 ng/ml	10 Tests	CE
FIA-CEA-001	CEA FIA Rapid Test Cassette	S/P	1~500 ng/ml	10 Tests	CE
FIA-CKMB-001	CK-MB FIA Rapid Test Cassette	WB/S/P	0.3-100 ng/ml	10 Tests	CE
FIA-CRP-001	CRP FIA Rapid Test Cassette	WB/S/P	0.5-200 mg/l	10 Tests	CE
FIA-CTI-001	cTnI FIA Rapid Test Cassette	WB/S/P	0.2-40 ng/ml	10 Tests	CE
FIA-CYSC-001	Cys-C FIA Rapid Test Cassette	WB/S/P	0.5-10 mg/l	10 Tests	CE
FIA-DDM-001	D-dimer FIA Rapid Test Cassette	WB/P	0.1~10 mg/l	10 Tests	CE
FIA-HBA-001	HbA1c FIA Rapid Test Cassette	WB	4-14.5%	10 Tests	CE
FIA-FABP-001	H-FABP FIA Rapid Test Cassette	WB/S/P	1-120 ng/ml	10 Tests	CE
FIA-MALB-001	mAlb FIA Rapid Test Cassette	Urine	5-300 mg/l	10 Tests	CE
FIA-MYO-001	Myoglobin FIA Rapid Test Cassette	WB/S/P	5-200 ng/ml	10 Tests	CE
FIA-BNP-001	NT-proBNP FIA Rapid Test Cassette	WB/S/P	18-35000 pg/ml	10 Tests	CE
FIA-PCT-001	PCT FIA Rapid Test Cassette	WB/S/P	0.1-50 ng/ml	10 Tests	CE
FIA-PSA-001	PSA FIA Rapid Test Cassette	S/P	2~100 ng/ml	10 Tests	Non-CE
FIA-STA-001	Strep A FIA Rapid Test Cassette	Swab	≥1E6	10 Tests	CE
FIA-TSH-001	TSH FIA Rapid Test Cassette	S/P	0.1-100 µIU/ml	10 Tests	CE
Fluorescence Immunoassay Analyser					
BR-FIA-1000	Fluorescence Immunoassay Analyser				CE



Biopanda Reagents
 Unit 14 Carrowreagh Business Park
 Carrowreagh Road, Belfast,
 BT16 1QQ, United Kingdom
 Email: info@biopanda.co.uk
 Web: www.biopanda.co.uk